



CO<sub>2</sub>

Dissolved carbon dioxide (DCO<sub>2</sub>) is a critical process parameter (CPP) in biopharma production processes according to PAT guidelines. By influencing other parameters such as extracellular and intracellular pH, it has an effect on different metabolic pathways which are involved in cell growth or in product formation and quality.

In the past, continuous in-line monitoring of DCO<sub>2</sub> has only been possible through electrochemical sensors that are based on the Severinghaus principle and measure the DCO<sub>2</sub> concentration indirectly. The result is significant maintenance effort and multiple sources of drift that must be compensated by time-consuming product calibration.

Now, Hamilton has introduced a completely new way to measure DCO<sub>2</sub>: The new in-line sensor CO<sub>2</sub>NTROL is a maintenance free, solid-state sensor that directly measures DCO<sub>2</sub> resulting in better measurement accuracy and lower cost of ownership.

# CO<sub>2</sub>NTROL **new**



CO<sub>2</sub>NTROL is the newest member to Hamilton's Arc Intelligent Sensor line. The Solid State Sensor directly measures DCO<sub>2</sub> and provides maintenance free, real-time, and in-line control of this new critical process parameter.

Unlike traditional sensors that are based on the electrochemical Severinghaus principle, CO<sub>2</sub>NTROL is a pure direct measurement in a solid state design: CO<sub>2</sub> molecules diffuse into a gas permeable membrane where the sensor measures the absorption of CO<sub>2</sub>-specific Mid-IR wavelengths. This absorption correlates to the partial pressure of CO<sub>2</sub> in the media.

CO<sub>2</sub>NTROL's hygienic design makes it compliant with requirements of biopharma applications. The sensor is EHEDG approved (EL Class I, test executed with Hamilton hygienic socket REF 242545) and is ready for GMP compliance. Embedded electronics convert the MIR CO<sub>2</sub> measurement into standard digital and analog signals that are easily integrated into your control strategy.

Arc Wi 2G Adapter BT (REF 243470) is required to output an analog 4-20 mA signal from the digital Modbus communication.



“**Did you know...**  
*Hamilton is the first and only supplier to bring the maintenance-free optical IR technology into a SIP/CIP compliant 12mm CO<sub>2</sub> sensor,*”

### Benefits

- ▶ Maintenance-free
- ▶ Simple calibration
- ▶ Hygienic design: SIP/CIP compatible, autoclavable
- ▶ Inverted installation possible
- ▶ Direct measurement of CO<sub>2</sub> – no ammonia interference

### Typical applications

- ▶ Biopharma Cell Cultures and Fermentations

### Specifications

Measurement Principle	Optical – CO <sub>2</sub> Absorption in Middle Infrared (MIR)
Measuring Range	5 to 1000 mbar or 0.5 to 100 %-Vol or 7.5 to 1500 mg/L (in liquid phase at 101.3 kPa and 25 °C)
Diameter	12 mm
Process Connection	PG 13.5
Wetted Parts	Stainless Steel 1.4435, EPDM (Ethylene propylene elastomer), FDA compliant silicone
Surface Quality	R <sub>a</sub> < 0.4 µm (N5)
Steam Sterilizable	Yes
Autoclavable	Yes
CIP	Yes
Operating temperature range	-10 to 60 °C

### Ordering Information

CO <sub>2</sub> NTROL	a-length	Arc
	120 mm	10087810-11
	160 mm	10087810-12
	225 mm*	10087810-13
	325 mm	10087810-14
	425 mm	10087810-15

\*CO<sub>2</sub>NTROL 225 have, in reality, a shaft length of 215 mm. This ensures optimal rinsing in replaceable armatures, such as Retractable.

### Accessories



**Calibration Station** Ref 243575

**Cables** see page 112

**Arc Accessories** see page 116

**Housings** see page 127

